## BT Scan

ZE/BTCAS25

ZE/BTCAS100

Rapid test for detection of betalactams  
and tetracyclines in milk

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## SCOPE AND TEST PRINCIPLE

**BT Scan** is a qualitative test designed for detection of betalactams and tetracyclines in milk. This test can be used in raw milk from cow, sheep and goat.

**BT Scan** is a competitive test which contains specific receptors and antibodies to detect betalactams ( $\beta$ ) and tetracyclines (T). If the sample is free from these antibiotics, the receptors will bind to their specific line of the strip when the liquid runs through the strip, and intense red lines will appear. When some antibiotic is present in the sample, it will bind the receptors which will not be able to bind to the specific line in the strip, inhibiting the appearance of the red color of the corresponding test line, or decreasing its intensity.

## KIT COMPONENTS

	ZE/BTCAS25	ZE/BTCAS100
TESTS	25	100
PIPETTE TIPS	25	100
PRODUCT CERTIFICATE	1	1
KIT INSTRUCTIONS	1	1

## ADDITIONAL MATERIAL (NOT PROVIDED)

- Micropipette 100  $\mu$ L
- Reader IRIS (reference ZE/IRIS)
- Smartphone

## SAFETY

A SAFETY DATA SHEET is available from your local distributor or ZEULAB. This kit should be used following good laboratory practices.

## STABILITY AND STORAGE

Store the test in dry conditions and under refrigeration (4-12°C) until its use. Check the expiry date on the package.

When the kit is used in trucks, it can be stored at room temperature (15-25°C) during the working day. Once finished the milk collection, the kit must be stored under refrigeration (4-12°C). It is recommended to use the kit in these conditions within a maximum period of 30 days.

## NOTES

- The test is very sensitive to moisture, do not open any bag until its use.
- Let the kit reach room temperature before its use.
- A new pipette tip must be used for each sample.
- The assay run automatically using IRIS reader. IRIS controls assay time and result interpretation.
- Please, contact ZEULAB for further information.

## TEST PROCEDURE

1. Access to the IRIS app in your smartphone.
2. Switch on the IRIS device and connect it to your smartphone app.
3. Wait until IRIS reaches its work temperature.
4. Select "Automatic assay".
5. Scan QR code from the product certificate.
6. Take out a cassette from the individual bag and place it in a preheated IRIS.
7. Add 100  $\mu$ L of milk sample into the circular cavity marked with an arrow on the cassette.
8. Press "Run test".
9. Read the results on your smartphone app. Assay takes 6 minutes.

**BT Scan** is an *in vitro* diagnostic kit. In analysis involving legal processes, the results should be re-evaluated with an official reference method. ZEULAB does not assume any legal responsibility. To the best of our knowledge, the information contained herein is accurate and complete. However, nothing herein shall be construed to imply any warranty or guarantee.

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