ECLIPSE 50

ZE/E50/96i
ZE/E50/288i

**Test for detection of inhibitory
substances in milk**

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SCOPE AND TEST PRINCIPLE

ECLIPSE is a qualitative test based on the inhibition of microbial growth, supplied in a microtiter plate format where each well contains agar medium spread with *Geobacillus stearothermophilus* spores and a pH indicator. When the plate is incubated at 65°C, spores germinate and cells grow producing acid and changing the agar pH. Variations of pH will produce changes of the agar colour from blue (purple) to yellowish. When milk samples contain inhibitors at higher concentrations than the detection limit, microorganisms will not grow and neither colour changes will be observed.

KIT COMPONENTS

Components	ZE/E50/96i	ZE/E50/288i
INDIVIDUAL TESTS	96	288
MICROTITER PLATE	1	3
ADHESIVE FILM	1	3
PRODUCT CERTIFICATE	Yes	Yes

ADDITIONAL MATERIAL (NOT PROVIDED)

- ✓ Micropipettes.
- ✓ Heater (FX incubator, ref: ZE/FX) or oven at 65°C.
- ✓ Negative control (sample without antibiotics)
- ✓ Positive control - freeze dried Penicillin G (ZE/PG5), oxytetracycline (ZE/OXITETRA) and sulfathiazole (ZE/SULFA).

SAFETY

Good laboratory practices are recommended when using this kit. A Material Safety Data Sheet (MSDS) is available from your local distributor or ZEULAB upon request.

STABILITY AND STORAGE

The kit components should be stored at 4-12 °C and protected from light. This kit is provided with a shelf life of approximately 9 months. See the kit expiry date on package.

NOTES

- ✓ A negative sample must be used to determine the optimal incubation of the assay in each run. A positive control sample is also recommended.
- ✓ A new pipette tip should be used for each sample.
- ✓ Sometimes milk can migrate to the bottom of the well. This will be identified as a white layer at the bottom of the well prior to washing. If you observe this phenomenon, you should read the results by looking all around the well sidewise.
- ✓ This test is extremely sensitive to antibiotics and other antibacterial substances, such as detergents and disinfectants. Any contamination with these substances should be prevented.
- ✓ Although natural inhibitors contained in milk do not interfere with the test results, samples from colostrum, milk from the end of the breeding period, and mastitis milk have high concentration of these inhibitors and can alter the results.
- ✓ Please, contact ZEULAB for the analysis of samples containing preservatives (i.e. Azidol).

TEST PROCEDURE

1. Cut the adhesive foil covering the wells and split the strips to be used by pressing up from the bottom. The foil covering the remaining wells should not be removed and the wells stored immediately at 4-12°C.
2. Remove the adhesive foil covering the wells/plate, and add 50 µl of sample into each well and seal carefully the wells with the adhesive film.
3. Incubate at 65°C until the negative control has turned to yellow (2h15'-2h45'). See the product certificate and use the incubation time shown as reference.
4. When the negative control has shifted to yellow as the color chart, incubation should be stopped. If it is not yet yellow, incubate the plate for longer (approx. 10-15min) until the negative control has changed. Test sensitivity may be affected increasing detection limits .
5. Remove the milk turning the plate and wash the wells with distilled water by filling them, empty them turning to the plate and tapping it against absorbent paper. Perform this washing step 2 or 3 times.
6. Results:
 - ✓ Visual reading: Turn the plate and read the results comparing with the negative control well.

- A blue-purple colour (positive) indicates the presence of antibiotics.
 - A green-blue colour (doubtful) indicates presence of antibiotics in a concentration close to the detection limit. Analysis should be repeated.
- ✓ **Photometric reading:** Read at 595nm (filter 1) and at 650nm (filter 2). The assay must be stopped when the difference of absorbance of negative control (AN 595nm - AN 650nm) is between 0.2 and 0.4 units. Mayor or equal values to those obtained for the samples used as negative control plus 0.2, will be positive.

$$\text{POSITIVE: AM 595nm - AM 650nm} \geq \text{AN 595nm - AN 650nm} + 0.2$$

AM: Sample absorbance

AN: Negative control absorbance

Note: The above procedure is only applicable when the difference of absorbance at 595nm - 650nm of the negative control sample is between 0.2 and 0.4.

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
Please contact ZEULAB or your local distributor for information on limit of detection of other species.

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SUMMARY PROCEDURE

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
Select the required number of tests for the analysis.

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
Incubación a 65°C

Wash and read:

- Visual reading from the bottom of the wells
- Photometrical reading at 595nm and 650nm

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-

-

ECLIPSE is an *in vitro* diagnostic kit for antibiotics screening in milk. In analysis implicating legal processes, the results should be revaluated with an official reference method. ZEULAB, S.L. does not assume any legal responsibility.

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