ECLIPSE FARM 3G

ZE/EF25

ZE/EF50

**Single tube tests for detection
of inhibitory substances in milk**

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info@zeulab.com • www.zeulab.com

SCOPE

ECLIPSE Farm^{3G} is a qualitative test designed for detection of antibiotics and inhibitors in milk. This test can be used in raw (with or without azidol), skimmed, heated or powder milk from cow, sheep, goat or buffalo. The analysis of raw cow milk can be performed in individual or commingled milk. For the other species, the analysis should be performed only with commingled milk.

PRINCIPLE

ECLIPSE Farm^{3G} is based on the inhibition of microbial growth. The tubes provided contain an specific culture medium spread with *Geobacillus stearothermophilus* spores and a pH indicator. When the tubes are incubated in absence of antibiotics, spores germinate and cells grow (negative). If samples contain inhibitors at higher concentrations than the detection limit, microorganisms will not grow (positive). The system has been optimized by combining the kit with e-Reader equipment, which integrates the processes of incubation, time control and interpretation of results.

KIT COMPONENTS

	ZE/EF25	ZE/EF50
Test tubes	25	50
Mini-pipettes	25	50
Adhesive film	1	2
Kit instructions	✓	✓
Product certificate	✓	✓

ADDITIONAL MATERIAL (NOT PROVIDED)

e-Reader (ref. ZE/ERT8) or FX incubator (ref. ZE/FX) or FS incubator (ref. ZE/FS).

Negative control sample (without antibiotics).

Positive control. Freeze dried Penicillin G (ref. ZE/PG5), Oxytetracycline (ref. ZE/OXITE-TRA), Sulfathiazole (ref. ZE/SULFA).

NOTES

- ✓ To achieve the best assay performance, it is recommended to analyse samples by using e-Reader, which incubates, stops the assay automatically and interprets the results (If it is not possible, use an FX or FS incubator, the procedure is described on the next page).
- ✓ Before starting, install and switch on e-Reader equipment following the manual provided.
- ✓ It is recommended to analyse a negative control sample of the same species and in the same conditions than the tested samples.
- ✓ A new pipette should be used for each sample. To use the provided pipettes: squeeze the bulb, put the stem in the sample and release the bulb to draw up liquid. Excess liquid is siphoned automatically into overflow tanks. Squeeze the bulb again to release 100 µl.

- ✓ This test is extremely sensitive to antibiotics and other antibacterial substances such as detergents and disinfectants. Any contamination with those substances should be prevented.
- ✓ Milk contains natural inhibitors which do not interfere with the result. However, these substances appear in high concentrations in colostrum, at the end of lactation and in mastitis cases, and they may alter the test results.
- ✓ The range of colours obtained using e-Reader is different to that seen with other incubation systems. Colours should not be compared with colour card shown on the last page.

TEST PROCEDURE (Flowchart Procedure on the last page)

Please pay attention to "Notes" in page 2.

1. Separate the tubes needed. Take care not to remove the foil on the remaining tubes to prevent the agar from drying out and stored immediately at 4-12°C.
2. Pierce or remove the foil covering each tube and add 100 µl of milk sample using the pipettes provided in the kit. (See "Notes" for pipettes use).
 - To detect antibiotics in cow milk: Seal carefully the test tubes with the adhesive film provided and place the sample tubes in e-Reader. Follow e-Reader instructions to begin the assay.
 - To detect antibiotics in goat or sheep milk, and samples containing azidiol*: Seal the test tubes with tape. Test tubes must be incubated at room temperature for 1 h. Afterwards, turn the tubes upside down to remove the remaining samples and wash the wells 2-3 times by filling them with distilled water. Remove the excess of water and seal carefully the test tubes with the adhesive film provided. Follow e-Reader instructions using Eclipse NC-SGA method.
- *Azidiol can be purchase from PanReac AppliChem (Product Code 176131).
3. e-Reader will stop automatically the assay and results will be shown on the screen.

Procedure without e-Reader

- Incubate the test tubes in a pre-warmed incubator at $65 \pm 1^\circ\text{C}$.
- Check the tubes at the time indicated on the product certificate. The incubation should be stopped when the negative control has turned yellow as shown in the colour card (from 2h15 to 3h). If the negative control sample is not yellow, incubate the tube again for further 10-15 min to get a clear yellow colour. When a negative control sample is not applied, the assay should be stopped after 3h of incubation. Please, be aware that sensitivity of the test can be affected, increasing the limit of detection in some antibiotics when stopping incubation at a fixed time of 3 hours.

- Results: (see colour chart)

- A yellow colour indicates the absence of antibiotics (negative).
- A blue-purple colour indicates the presence of antibiotics (positive).
- A green-blue colour indicates presence of antibiotics in a concentration close to the detection limit (doubtful). Analysis should be repeated.

STABILITY AND STORAGE

The kit components should be stored at 4-12°C and protected from light. The kit is provided with a minimum shelf life of 9 months. See the kit expiring date on the package.

SAFETY

Good laboratory practices are recommended when using this kit. A Material Safety Sheet (MSDS) is available from your local distributor or ZEULAB upon request.

ECLIPSE is an in vitro diagnostic kit for antibiotics screening in milk. In analysis implicating legal processes, the results should be reevaluated with an official reference method. ZEULAB does not assume any legal responsibility.

To the best of our knowledge, the information contained herein is accurate and complete. However, nothing herein shall be construed to imply any warranty or guarantee.

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